

Anti-Znf143 Rabbit Polyclonal Antibody (Wb, Elisa) - P52747

Item Number: CM0019858

Introduction

Rabbit polyclonal antibody targeting human ZNF143 (SPH-binding factor/hStaf) for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for research on transcriptional regulation and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	ZNF143 Rabbit pAb
Remarks/Alias	SBF; STAF; pHZ-1; ZNF143
Species	Human
GeneID	7702
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 437-638 of human ZNF143 (NP_003433.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52747
Protein Weight	69kDa
Shipping	Ice bag