

Anti-Znf16 (Hzf1/Kox9) Polyclonal Antibody For Wb And Elisa - P17020

Item Number: CM0021011

Introduction

Rabbit polyclonal antibody targeting human ZNF16 (HZF1/KOX9), a zinc-finger transcription factor involved in cell cycle and apoptosis. Suitable for WB and ELISA. Reacts with human species. 76kDa.

[Learn More](#)

Parameter	Value
Product Name	ZNF16 Rabbit pAb
Remarks/Alias	HZF1, KOX9, ZNF16
Species	Human
Gene ID	7564 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human ZNF16 (NP_008889.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17020
Protein Weight	76kDa
Shipping	Ice bag