

Anti-Znf177 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q13360

Item Number: CM0014007

Introduction

Rabbit polyclonal antibody targeting human ZNF177, validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for transcriptional regulation studies in research settings.

[Learn More](#)

Parameter	Value
Product Name	ZNF177 Rabbit pAb
Remarks/Alias	PIGX; ZNF177
Species	Human
Gene ID (Human)	7730
Gene ID	7730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-110 of human ZNF177 (NP_003442.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13360
Protein Weight	55kDa
Shipping	Ice bag