

Anti-Znf182 Polyclonal Antibody For Wb And Elisa - P17025

Item Number: CM0023898

Introduction

High-quality Anti-ZNF182 rabbit polyclonal antibody for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat species. Targets the zinc finger protein 182 transcriptional regulator.

[Learn More](#)

Parameter	Value
Product Name	ZNF182 Rabbit pAb
Remarks/Alias	KOX14; ZNF21; HHZ150; Zfp182; ZNF182
Species	Human
Gene ID (Human)	7569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human ZNF182 (NP_001171570.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17025
Protein Weight	74kDa
Shipping	Ice bag