

Anti-Znf195 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O14628

Item Number: CM0026020

Introduction

Rabbit polyclonal antibody targeting ZNF195 (SWISS O14628), validated for WB, IF/ICC, ELISA. Detects human, mouse, rat. Ideal for studying transcription regulation and zinc finger protein biology.

[Learn More](#)

Parameter	Value
Product Name	ZNF195 Rabbit pAb
Remarks/Alias	HRF1; ZNFP104; ZNF195
Species	Human
Gene ID (Human)	7748
Gene ID	7748
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ZNF195 (NP_001123991.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14628
Protein Weight	72kDa
Shipping	Ice bag