

Anti-Znf20 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17024

Item Number: CM0030654

Introduction

ZNF20 Rabbit pAb (human ZNF20/KOX13, gene ID 7568) for Western blot, IHC-P, ELISA. Validated in human, mouse, rat. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF20 Rabbit pAb
Remarks/Alias	KOX13; ZNF20
Species	Human
Gene ID	7568 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ZNF20 (NP_066966.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17024
Protein Weight	62kDa
Shipping	Ice bag