

Anti-Znf20 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P17024

Item Number: CM0007208

Introduction

Recombinant rabbit monoclonal antibody targeting human zinc finger protein 20 (ZNF20/KOX13), validated for WB, IF-P, IHC-P, and ELISA applications. Detects mouse and rat orthologs. Useful for transcriptional regulation studies.

[Learn More](#)

Parameter	Specification
Product Name	ZNF20 Rabbit PolymAb®
Remarks/Alias	KOX13
Species	Human
Gene ID (human)	7568
Gene ID	7568
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 51-118 of human ZNF20 (NP_066966.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P17024
Protein Weight	62kDa
Shipping	Ice bag