

Anti-Znf207 (Bugz, Hbugz) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43670

Item Number: CM0023129

Introduction

Rabbit polyclonal antibody against human ZNF207 (BuGZ, hBuGZ), validated for WB, IHC-P, ELISA; cross-reacts with mouse, rat. For spindle assembly and chromosome alignment research.

[Learn More](#)

Parameter	Value
Product Name	ZNF207 Rabbit pAb
Remarks/Alias	BuGZ; hBuGZ; ZNF207
Species	Human
Gene ID	7756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human ZNF207 (NP_003448.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43670
Protein Weight	51kDa
Shipping	Ice bag