

Anti-Znf213 Polyclonal Antibody For Wb, Ihc-P, Elisa - O14771

Item Number: CM0021792

Introduction

Rabbit polyclonal antibody targeting human ZNF213 (CR53/ZKSCAN21), validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for transcriptional regulation research and nuclear protein studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF213 Rabbit pAb
Remarks/Alias	CR53; ZSCAN53; ZKSCAN21; ZNF213
Species	Human
GeneID (Human)	7760
GeneID	7760
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 132-459 of human ZNF213 (NP_004211.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14771
Protein Weight	51kDa
Shipping	Ice bag