

Anti-Znf22 Rabbit Polyclonal Antibody For Wb, Elisa - P17026

Item Number: CM0014170

Introduction

ZNF22 rabbit polyclonal antibody validated for WB and ELISA. Detects human ZNF22 and cross-reacts with mouse. Ideal for transcriptional regulation studies and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	ZNF22 Rabbit pAb
Remarks/Alias	KOX15; HKR-T1; ZNF422; Zfp422; ZNF22
Species	Human
Gene ID (Human)	7570
Gene ID	7570
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human ZNF22 (NP_008894.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P17026
Protein Weight	26kDa
Shipping	Ice bag