

Anti-Znf23 Rabbit Polyclonal Antibody For Wb, Elisa - P17027

Item Number: CM0014443

Introduction

High-quality rabbit polyclonal antibody targeting human ZNF23 (KOX16, ZNF359, ZNF612), validated for Western blot and ELISA. Shows cross-reactivity with mouse. Suitable for studies in transcriptional regulation and embryonic development.

[Learn More](#)

Parameter	Value
Product Name	ZNF23 Rabbit pAb
Remarks/Alias	KOX16; ZNF359; ZNF612; Zfp612; ZNF23
Species	Human
GeneID (Human)	7571
GeneID	7571
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ZNF23 (NP_666016.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17027
Protein Weight	73kDa
Shipping	Ice bag