

Anti-Znf24 Polyclonal Antibody For Wb, Elisa - P17028

Item Number: CM0028196



Introduction

High-quality rabbit anti-ZNF24 polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat ZNF24 (KOX17/RSG-A). Suitable for transcription factor and myelination studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF24 Rabbit pAb
Remarks/Alias	KOX17; RSG-A; ZNF191; ZSCAN3; Zfp191; ZNF24
Species	Human
GeneID (Human)	7572
GeneID	7572
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 169-368 of human ZNF24 (NP_008896.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P17028
Protein Weight	42kDa
Shipping	Ice bag