

Anti-Znf248 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ndw4

Item Number: CM0012977

Introduction

Anti-ZNF248 rabbit polyclonal antibody validated for Western blot and ELISA. Targets human zinc finger protein 248, a putative transcriptional regulator. Suitable for transcription research and protein detection.

[Learn More](#)

Parameter	Value
Product Name	ZNF248 Rabbit pAb
Remarks/Alias	bA162G10.3; ZNF248
Species	Human
GeneID	57209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human ZNF248 (NP_066383.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NDW4
Protein Weight	67kDa
Shipping	Ice bag