

Anti-Znf259 Monoclonal Antibody For Wb, Ihc-P, Elisa - 075312

Item Number: CM0003059

Introduction

Rabbit monoclonal antibody against ZNF259 (ZPR1), validated for WB, IHC-P, ELISA. Reacts with human and mouse samples. Ideal for studying cell cycle regulation and SMN1 localization.

[Learn More](#)

Parameter	Value
Product Name	ZNF259 Rabbit mAb
Remarks/Alias	GKAF; ZNF259
Species	Human
Gene ID (Human)	8882
Gene ID	8882
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 360-459 of human ZNF259 (075312).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	075312
Protein Weight	51kDa
Shipping	Ice bag