

Anti-Znf263 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14978

Item Number: CM0022469

Introduction

High-quality rabbit polyclonal antibody specific for ZNF263, a KRAB/SCAN zinc finger transcription factor that represses gene expression via chromatin modification. Validated for Western blot, immunohistochemistry, and ELISA in human, mouse, and rat samples. Suitable for transcription regulation and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	ZNF263 Rabbit pAb
Remarks/Alias	FPM315; ZSCAN44; ZKSCAN12; ZNF263
Species	Human
GeneID (Human)	10127
GeneID	10127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-240 of human ZNF263 (NP_005732.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14978
Protein Weight	77kDa
Shipping	Ice bag