

Anti-Znf264 Polyclonal Antibody For Wb, Elisa - 043296

Item Number: CM0021162

Introduction

High-quality rabbit polyclonal antibody targeting human ZNF264 (O43296). Validated for Western blot (WB) and ELISA. Ideal for transcriptional regulation, DNA binding, and nuclear protein research. Reliable IVD raw material sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	ZNF264 Rabbit pAb
Species	Human
Gene ID (Human)	9422
Gene ID	9422
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human ZNF264. (NP_003408.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43296
Protein Weight	71kDa
Shipping	Ice bag