

Anti-Znf268 (Hzf3) Rabbit Polyclonal Antibody For Wb, Elisa - Q14587

Item Number: CM0012962



Introduction

Rabbit polyclonal antibody targeting ZNF268 (HZF3), a transcriptional repressor involved in cancer and antiviral signaling. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for transcription regulation and oncology research.

[Learn More](#)

Parameter	Value
Product Name	ZNF268 Rabbit pAb
Remarks/Alias	HZF3; ZNF268
Species	Human
GeneID (human)	10795
GeneID	10795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ZNF268 (NP_001159353.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14587
Protein Weight	108kDa
Shipping	Ice bag