

Anti-Znf273 Rabbit Polyclonal Antibody For Wb, Elisa - Q14593

Item Number: CM0015795

Introduction

Rabbit polyclonal antibody to human ZNF273 (HZF9), validated for WB and ELISA. A zinc finger transcription factor localized to the nucleus, suitable for transcriptional regulation research.

[Learn More](#)

Parameter	Specification
Product Name	ZNF273 Rabbit pAb
Remarks/Alias	HZF9; ZNF273
Species	Human
GeneID (Human)	10793
GeneID	10793
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 210-260 of human ZNF273 (NP_066971.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14593
Protein Weight	65kDa
Shipping	Ice bag