

Anti-Znf296 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wuu4

Item Number: CM0014609

Introduction

ZNF296 rabbit polyclonal antibody validated for WB and ELISA, reactive with human and mouse. Ideal for studying transcriptional regulation, zinc finger proteins, and KLF4 interactions. Research-grade IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	ZNF296 Rabbit pAb
Remarks/Alias	ZFP296; ZNF342; ZNF296
Species	Human
Gene ID (Human)	162979
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ZNF296 (NP_660331.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WUU4
Protein Weight	51kDa
Shipping	Ice bag