

Anti-Znf302 Polyclonal Antibody For Wb, Elisa - Q9Nr11

Item Number: CM0020400

Introduction

Rabbit polyclonal antibody targeting human ZNF302 (zinc finger protein 302), suitable for Western blot and ELISA with cross-reactivity to rat. Ideal for transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF302 Rabbit pAb
Remarks/Alias	HSD16; MST154; ZNF327; MSTP154; ZNF135L; ZNF140L; ZNF302
Species	Human
Gene ID (Human)	55900
Gene ID	55900
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-200 of human ZNF302 (NP_060913.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9NR11
Protein Weight	55kDa
Shipping	Ice bag