

Anti-Znf324 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - 075467

Item Number: CM0020726

Introduction

ZNF324 Rabbit pAb for human ZNF324, validated in IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for transcriptional regulation research and cell proliferation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF324 Rabbit pAb
Remarks/Alias	ZNF324; ZF5128; ZNF324A; zinc finger protein 324
Species	Human
GeneID (Human)	25799
GeneID	25799
Immunogen	Recombinant protein corresponding to a sequence within amino acids 271-553 of human ZNF324 (NP_055162.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75467
Protein Weight	61kDa
Shipping	Ice bag