

Anti-Znf333 Rabbit Polyclonal Antibody For Wb, Elisa - Q96JL9

Item Number: CM0017766

Introduction

Rabbit polyclonal antibody targeting human ZNF333. Validated for Western blot and ELISA. Cross-reacts with mouse and rat ZNF333. Suitable for transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF333 Rabbit pAb
Remarks/Alias	ZNF333
Species	Human
Gene ID	84449
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-303 of human ZNF333 (NP_001287841.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96JL9
Protein Weight	76 kDa
Shipping	Ice bag