

Anti-Znf346 Rabbit Polyclonal Antibody For Wb And Elisa - Q9UI40

Item Number: CM0025993



Introduction

Rabbit polyclonal antibody targeting human ZNF346 (JAZ zinc finger protein). Suitable for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Recombinant protein immunogen. Ideal for RNA-binding and nucleolar function research.

[Learn More](#)

Parameter	Value
Product Name	ZNF346 Rabbit pAb
Remarks/Alias	JAZ; Zfp346; ZNF346
Species	Human
GeneID (Human)	23567
GeneID	23567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-294 of human ZNF346 (NP_036411.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL40
Protein Weight	33kDa
Shipping	Ice bag