

Anti-Znf350 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Gzx5

Item Number: CM0029952

Introduction

Rabbit polyclonal antibody against human ZNF350 (ZBRK1, ZFQR) validated for WB and ELISA. Detects human and mouse samples. Suitable for studies on transcriptional repression and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	ZNF350 Rabbit pAb
Remarks/Alias	ZFQR; ZBRK1; ZNF350
Species	Human
Gene ID (Human)	59348
Gene ID	59348
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 285-532 of human ZNF350 (NP_067645.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9GZX5
Protein Weight	60kDa
Shipping	Ice bag