

Anti-Znf408 Rabbit Pab For Wb, Elisa - Q9H9D4

Item Number: CM0027463

Introduction

High-quality rabbit polyclonal antibody against human ZNF408, validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for transcriptional regulation and retinal disease research.

[Learn More](#)

Parameter	Value
Product Name	ZNF408 Rabbit pAb
Remarks/Alias	EVR6; RP72; ZNF408
Species	Human
GeneID	79797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 471-720 of human ZNF408 (NP_079017.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9D4
Protein Weight	78kDa
Shipping	Ice bag