

## Anti-Znf416 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q9Bwm5

Item Number: CM0017953

### Introduction

CamelBio provides anti-ZNF416 rabbit polyclonal antibody suitable for IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat ZNF416, a mechanosensitive transcriptional activator in fibroblast biology research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | ZNF416 Rabbit pAb                                                                                                   |
| Remarks/Alias    | ZNF416                                                                                                              |
| Species          | Human                                                                                                               |
| GeneID           | 55659                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 495-594 of human ZNF416 (NP_060349.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | IHC-P, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q9BWM5                                                                                                              |
| Protein Weight   | 67kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |