

Anti-Znf433 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q8N7K0

Item Number: CM0015573

Introduction

Anti-ZNF433 rabbit polyclonal antibody validated for IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat ZNF433, a putative zinc-finger transcription regulator. Suitable for gene regulation and nuclear protein studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF433 Rabbit pAb
Remarks/Alias	ZNF433
Species	Human
GeneID (Human)	163059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 564-673 of human ZNF433 (NP_001073880.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N7K0
Protein Weight	77kDa
Shipping	Ice bag