

Anti-Znf451 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4E5

Item Number: CM0013941



Introduction

Rabbit polyclonal antibody against human ZNF451 (E3 SUMO-protein ligase, also known as COASTER/ZATT). Validated for Western blot and ELISA in human, mouse, and rat. Ideal for studying SUMOylation, DNA damage response, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	ZNF451 Rabbit pAb
Remarks/Alias	ZATT; COASTER; dj417l1.1; ZNF451
Species	Human
GeneID (Human)	26036
GeneID	26036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-900 of human ZNF451 (NP_001026794.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4E5
Protein Weight	121kDa
Shipping	Ice bag