

Anti-Znf460 Polyclonal Antibody For Wb And Elisa - Q14592

Item Number: CM0022015

Introduction

ZNF460 rabbit polyclonal antibody (alias HZF8, ZNF272) validated for Western blot and ELISA. Detects human, mouse, and rat ZNF460. Suitable for nuclear transcription factor studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF460 Rabbit pAb
Remarks/Alias	HZF8; ZNF272; ZNF460
Species	Human
GeneID (Human)	10794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human ZNF460 (NP_006626.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14592
Protein Weight	64kDa
Shipping	Ice bag