

Anti-Znf471 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bx82

Item Number: CM0024354

Introduction

ZNF471 rabbit polyclonal antibody validated for WB and ELISA. Specific to human ZNF471, a zinc finger transcription factor. Ideal for epigenetic and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF471 Rabbit pAb
Remarks/Alias	ERP1; Z1971; Zfp78; ZNF471
Species	Human
GeneID (human)	57573
GeneID	57573
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-221 of human ZNF471 (NP_065864.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BX82
Protein Weight	73kDa
Shipping	Ice bag