

Anti-Znf530 Rabbit Polyclonal Antibody For Wb, Elisa - Q6P9A1

Item Number: CM0024661

Introduction

Rabbit polyclonal antibody against human ZNF530, a likely transcriptional regulator. Validated for Western blot and ELISA; cross-reacts with human, mouse, and rat. Suitable for gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	ZNF530 Rabbit pAb
Remarks/Alias	ZNF530
Species	Human
Gene ID (Human)	348327
Gene ID	348327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ZNF530 (NP_065931.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P9A1
Protein Weight	69kDa
Shipping	Ice bag