

Anti-Znf562 Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q6V9R5

Item Number: CM0015605

Introduction

Rabbit polyclonal antibody targeting ZNF562, a zinc finger protein involved in transcriptional regulation. Verified for IHC-P and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for gene expression and nuclear localization studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF562 Rabbit pAb
Remarks/Alias	ZNF562
Species	Human
GeneID (Human)	54811
GeneID	54811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ZNF562 (NP_001123503.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6V9R5
Protein Weight	49kDa
Shipping	Ice bag