

Anti-Znf563 Polyclonal Antibody For Wb And Elisa - Q8Ta94

Item Number: CM0019279

Introduction

ZNF563 rabbit polyclonal antibody for Western blot (WB) and ELISA applications. Validated in human, mouse, and rat species. Target: zinc finger protein 563, a predicted transcription factor. Ideal for transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF563 Rabbit pAb
Remarks/Alias	ZNF563
Species	Human
Gene ID (Human)	147837
Gene ID	147837
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ZNF563 (NP_660319.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TA94
Protein Weight	55kDa
Shipping	Ice bag