

## Anti-Znf566 Rabbit Polyclonal Antibody For Wb And Elisa - Q969W8

Item Number: CM0020027

### Introduction

Rabbit polyclonal antibody against human ZNF566, validated for Western blot and ELISA. Designed for transcriptional regulation studies; detects a 49 kDa nuclear zinc finger protein. Ideal for research applications.

[Learn More](#)

| Parameter        | Value                                                                                                      |
|------------------|------------------------------------------------------------------------------------------------------------|
| Product Name     | ZNF566 Rabbit pAb                                                                                          |
| Remarks/Alias    | ZNF566; zinc finger protein 566; 66                                                                        |
| Species          | Human                                                                                                      |
| GeneID           | 84924                                                                                                      |
| Immunogen        | Recombinant protein corresponding to a sequence within amino acids 50-154 of human ZNF566 (NP_001138816.1) |
| Source           | Rabbit                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                      |
| Application      | WB, ELISA                                                                                                  |
| Cross Reactivity | Human                                                                                                      |
| SWISS            | Q969W8                                                                                                     |
| Protein Weight   | 49kDa                                                                                                      |
| Shipping         | Ice bag                                                                                                    |