

Anti-Znf574 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q6Zn55

Item Number: CM0014495

Introduction

Rabbit polyclonal antibody against ZNF574, a zinc finger protein involved in ribosome assembly surveillance. Validated for IHC-P, IF/ICC, ELISA in human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	ZNF574 Rabbit pAb
Remarks/Alias	FP972; ZNF574
Species	Human
GeneID (human)	64763
GeneID	64763
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-310 of human ZNF574 (NP_073589.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZN55
Protein Weight	99kDa
Shipping	Ice bag