

Anti-Znf582 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ng8

Item Number: CM0015036

Introduction

Anti-ZNF582 Rabbit Polyclonal Antibody for WB and ELISA, with mouse cross-reactivity. Targets ~60kDa nuclear zinc-finger protein likely involved in transcriptional regulation, suitable for gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF582 Rabbit pAb
Remarks/Alias	ZNF582
Species	Human
GeneID (Human)	147948
GeneID	147948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-180 of human ZNF582 (NP_653291.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96NG8
Protein Weight	60kDa
Shipping	Ice bag