

Anti-Znf587 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Sq5

Item Number: CM0023444

Introduction

Rabbit polyclonal antibody against human ZNF587, a zinc finger transcriptional regulator. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	ZNF587 Rabbit pAb
Remarks/Alias	ZF6; ZNF587
Species	Human
Gene ID	84914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-240 of human ZNF587 (NP_116217.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SQ5
Protein Weight	66kDa
Shipping	Ice bag