

Anti-Znf598 Rabbit Polyclonal Antibody For Wb And Elisa - Q86Uk7

Item Number: CM0018977

Introduction

Rabbit polyclonal antibody targeting ZNF598, an E3 ubiquitin ligase in ribosome quality control. Validated for Western blot and ELISA with human samples. Ideal for translation regulation and viral infection studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF598 Rabbit pAb
Remarks/Alias	HEL2; ZNF598
Species	Human
Gene ID (Human)	90850
Gene ID	90850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-420 of human ZNF598 (NP_835461.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86UK7
Protein Weight	99kDa
Shipping	Ice bag