

Anti-Znf620 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Zng0

Item Number: CM0019701

Introduction

Rabbit polyclonal antibody targeting human ZNF620, validated for WB and ELISA. Cross-reacts with mouse. Predicted molecular weight ~49 kDa. Suitable for nuclear transcription factor studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF620 Rabbit pAb
Remarks/Alias	ZNF620
Species	Human
GeneID (Human)	253639
GeneID	253639
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human ZNF620 (NP_001243096.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q6ZNG0
Protein Weight	49kDa
Shipping	Ice bag