

Anti-Znf621 Polyclonal Antibody For Wb, Elisa - Q6Zss3

Item Number: CM0013816

Introduction

Rabbit polyclonal antibody targeting human ZNF621, a zinc finger transcription regulator. Validated for Western Blot and ELISA applications. Suitable for gene expression and protein interaction studies. - Q6ZSS3

[Learn More](#)

Parameter	Value
Product Name	ZNF621 Rabbit pAb
Remarks/Alias	ZNF621
Species	Human
GeneID (Human)	285268
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-150 of human ZNF621 (NP_940886.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6ZSS3
Protein Weight	49kDa
Shipping	Ice bag