

Anti-Znf645 Polyclonal Antibody For Wb - Q8N7E2

Item Number: CM0019273

Introduction

Rabbit polyclonal antibody against human ZNF645 (E3 ubiquitin-protein ligase CBLL2), validated for WB and ELISA. Detects mouse and rat proteins. Essential for ubiquitin ligase and testis biology studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF645 Rabbit pAb
Remarks/Alias	CT138, HAKAIL, ZNF645
Species	Human
GeneID (Human)	158506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 346-425 of human ZNF645 (NP_689790.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N7E2
Protein Weight	49kDa
Shipping	Ice bag