

Anti-Znf652 Polyclonal Antibody For Wb, Elisa - Q9Y2D9

Item Number: CM0022935

Introduction

Rabbit polyclonal antibody against human ZNF652, suitable for WB and ELISA. Cross-reacts with mouse and rat. Recognizes a ~70 kDa protein. Ideal for transcription repression studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF652 Rabbit pAb
Remarks/Alias	ZNF652
Species	Human
Gene ID (Human)	22834
Gene ID	22834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human ZNF652 (NP_055712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2D9
Protein Weight	70kDa
Shipping	Ice bag