

Anti-Znf652 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y2D9

Item Number: CM0023648

Introduction

CamelBio supplies this rabbit anti-ZNF652 polyclonal antibody, validated for WB and ELISA, targeting the transcription repressor ZNF652 across human, mouse, and rat samples. Ideal for gene regulation and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	ZNF652 Rabbit pAb
Species	Human
Gene ID (Human)	22834
Gene ID	22834
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human ZNF652 (NP_055712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2D9
Protein Weight	70kDa
Shipping	Ice bag