

Anti-Znf677 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Xu0

Item Number: CM0024315

Introduction

Rabbit polyclonal antibody targeting human ZNF677 (68 kDa) for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for studying transcriptional regulation and DNA-binding nuclear proteins.

[Learn More](#)

Parameter	Value
Product Name	ZNF677 Rabbit pAb
Remarks/Alias	ZNF677
Species	Human
GeneID (Human)	342926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human ZNF677 (NP_872415.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86XU0
Protein Weight	68kDa
Shipping	Ice bag