

Anti-Znf683 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Iz20

Item Number: CM0027090

Introduction

Rabbit polyclonal antibody against human ZNF683 (Hobit) for Western blot and ELISA. Predicted molecular weight 57 kDa. Suitable for studying tissue-resident lymphocyte biology and immune cell retention.

[Learn More](#)

Parameter	Value
Product Name	ZNF683 Rabbit pAb
Remarks/Alias	Hobit, ZNF683
Species	Human
Gene ID (Human)	257101
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 420-524 of human ZNF683 (NP_001294854.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IZ20
Protein Weight	57kDa
Shipping	Ice bag