

Anti-Znf688 Polyclonal Antibody For Wb, Elisa - P0C7X2

Item Number: CM0016151

Introduction

Rabbit polyclonal antibody targeting ZNF688, a zinc finger protein potentially involved in transcriptional regulation, validated for Western blot (WB) and ELISA in human, mouse, and rat samples. Ideal for gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF688 Rabbit pAb
Remarks/Alias	ZNF688
Species	Human
GeneID (Human)	146542
GeneID	146542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human ZNF688 (NP_660314.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0C7X2
Protein Weight	31kDa
Shipping	Ice bag