

Anti-Znf711 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y462

Item Number: CM0020804



Introduction

ZNF711 rabbit polyclonal antibody validated for WB and ELISA. Targets human ZNF711 (ZNF6, CMPX1), a transcription regulator essential for brain development. Ideal for neuroscience research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	ZNF711 Rabbit pAb
Remarks/Alias	ZNF4; ZNF5; ZNF6; CMPX1; MRX65; MRX97; XLID97; Zfp711; dj75N13.1; ZNF711
Species	Human
GeneID (Human)	7552
GeneID	7552
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ZNF711 (NP_001317503.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y462
Protein Weight	86kDa
Shipping	Ice bag