

## Anti-Znf75A Rabbit Polyclonal Antibody For Wb, Elisa - Q96N20

Item Number: CM0020364

### Introduction

ZNF75A rabbit polyclonal antibody validated for WB and ELISA. Recognizes human ZNF75A, a transcriptional regulator. For gene regulation research.

[Learn More](#)

| Parameter        | Value                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------------|
| Product Name     | ZNF75A Rabbit pAb                                                                                        |
| Remarks/Alias    | ZNF75A                                                                                                   |
| Species          | Human                                                                                                    |
| GeneID (human)   | 7627                                                                                                     |
| GeneID           | 7627                                                                                                     |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 200-296 of human ZNF75A (NP_694573.1) |
| Source           | Rabbit                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                    |
| Application      | WB, ELISA                                                                                                |
| Cross Reactivity | Human                                                                                                    |
| SWISS            | Q96N20                                                                                                   |
| Protein Weight   | 35kDa                                                                                                    |
| Shipping         | Ice bag                                                                                                  |