

Anti-Znf76 Polyclonal Antibody For Wb And Elisa - P36508

Item Number: CM0024059

Introduction

Rabbit polyclonal antibody targeting ZNF76, a transcription factor regulating RNA pol II/III. Validated for WB, ELISA; cross-reacts with mouse and rat. Ideal for gene expression and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	ZNF76 Rabbit pAb
Remarks/Alias	ZNF76; D6S229E; ZNF523; Zfp523; zinc finger protein 76
Species	Human
Gene ID	7629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 463-570 of human ZNF76 (NP_003418.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P36508
Protein Weight	56kDa/61kDa
Shipping	Ice bag