

Anti-Znf76 Rabbit Polyclonal Antibody For Wb, Elisa - P36508

Item Number: CM0022840

Introduction

High-quality rabbit polyclonal antibody against human ZNF76 transcription factor for Western blot and ELISA applications. Validated in human, mouse, and rat samples. Research-grade IVD raw material.

[Learn More](#)

Parameter	Specification
Product Name	ZNF76 Rabbit pAb
Remarks/Alias	ZNF523; Zfp523; D6S229E; ZNF76
Species	Human
GeneID (Human)	7629
GeneID	7629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-170 of human ZNF76 (NP_003418.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36508
Protein Weight	62kDa
Shipping	Ice bag